

LIPO-SUCTION CANNULAS

07



PAGE 001 to 009



- EN Lipo-Suction Cannulas
- DE Fettabsaugung Cannulas
- ES Cánulas del Lipo-Suction

SERVING HEALTHCARE WORLDWIDE

www.hilbro.net



Lipo-Suction Cannulas EN

Fettabsaugung Cannulas DE

Cánulas del Lipo-Suction ES

07

LIPO-SUCTION CANNULAS

Hilbro



One central hole,
two lateral holes



Two lateral holes,
one central hole



Three holes, diagonal



	Ø	L		Ø	L		Ø	L
07.0030.20	3.0mm	20.0cm/8"	07.0052.20	3.0mm	20.0cm/8"	07.0074.20	3.0mm	20.0cm/8"
07.0030.30	3.0mm	30cm/12"	07.0052.30	3.0mm	30cm/12"	07.0074.30	3.0mm	30cm/12"
07.0034.30	4.0mm	30cm/12"	07.0056.30	4.0mm	30cm/12"	07.0078.30	4.0mm	30cm/12"
07.0038.30	5.0mm	30cm/12"	07.0060.30	5.0mm	30cm/12"	07.0078.40	4.0mm	40cm/15¾"
07.0038.40	5.0mm	40cm/15¾"	07.0064.30	6.0mm	30cm/12"	07.0082.30	5.0mm	30cm/12"
07.0042.30	6.0mm	30cm/12"				07.0086.30	6.0mm	30cm/12"

L = Working Length

07.0086.30

07-001

www.hilbro.net

07

EN Lipo-Suction Cannulas

DE Fettabsaugung Cannulas

ES Cánulas del Lipo-Suction

LIPO-SUCTION CANNULAS



One central hole,
two lateral holes

Tip with hole,
one central hole

10 x 6mm

	Ø	L		Ø	L		Ø	L
07.0096.30	4.0mm	30cm/12"	07.0110.20	3.0mm	20.0cm/8"	07.0120.30	6.0mm	30.0cm/12"
07.0100.30	5.0mm	30cm/12"	07.0110.30	3.0mm	30cm/12"	07.0120.35	6.0mm	35cm/13¾"


07.0120.35

07-002

www.hilbro.net

L = Working Length

LIPO-SUCTION CANNULAS



	Ø	L		Ø	L		Ø	L
07.0130.05	2.0mm	5.0 cm/2"	07.0172.15	3.0mm	15.0cm/6"	07.0186.30	3.0mm	30cm/12"
07.0134.05	3.0mm	5.0 cm/2"	07.0172.25	3.0mm	25cm/10"			
07.0138.05	4.0mm	5.0 cm/2"	07.0176.15	4.0mm	15.0cm/6"			
07.0142.15	2.0mm	15.0cm/6"	07.0176.25	4.0mm	25cm/10"			
07.0146.15	3.0mm	15.0cm/6"						
07.0146.25	3.0mm	25cm/10"						
07.0150.25	4.0mm	25cm/10"						
07.0154.30	5.0mm	30cm/12"						
07.0158.30	3.0mm	30cm/12"						
07.0162.30	4.0mm	30cm/12"						

L = Working Length

07.0186.30

07-003

www.hilbro.net

07

EN Lipo-Suction Cannulas


DE Fettabsaugung Cannulas

ES Cánulas del Lipo-Suction


LIPO-SUCTION CANNULAS

For fine facial contourings

For fine facial & neck contourings



One central hole



8 holes


Handle conn. Ø	(L)	Luer-Lock	Handle conn. Ø	Ø Hole	(L)	Luer-Lock
07.0196.05	1.7mm	05cm/2"	07.2196.05	1.5mm	0.5mm	10cm/4"
07.0200.05	2.0mm	05cm/2"	07.2200.05	2.0mm	1.0mm	10cm/4"
07.0204.05	3.0mm	05cm/2"	07.2204.05	2.0mm	1.0mm	15cm/6"
07.0208.15	1.7mm	15cm/6"	07.2208.15			
07.0212.15	2.0mm	15cm/6"	07.2212.15			
07.0216.15	3.0mm	15cm/6"	07.2216.15			

07.2230.15

07-004

www.hilbro.net

L = Working Length

Lipo-Suction Cannulas **EN**Fettabsaugung Cannulas **DE**Cánulas del Lipo-Suction **ES****07****LIPO-SUCTION CANNULAS**

Handle conn.	Ø	(L)	Hole (Ø Hole)	Luer-Lock	Handle conn.	Ø	(L)	Hole (Ø Hole)	Luer-Lock
07.0240.15	2.0mm	15.0cm/6"	18 1.0mm	07.2240.15	07.0278.20	3.0mm	20 cm/8"	22 1.5mm	07.2278.20
07.0241.15	2.5mm	15.0cm/6"	18 1.2mm		07.0282.20	4.0mm	20 cm/8"	22 2.0mm	07.2282.20
07.0244.15	3.0mm	15.0cm/6"	18 1.5mm	07.2244.15	07.0286.30	3.0mm	30cm/12"	30 1.5mm	07.2286.30
07.0249.20	2.5mm	20.0cm/8"	22 1.2mm		07.0290.30	4.0mm	30cm/12"	30 2.0mm	07.2290.30
07.0252.20	3.0mm	20.0cm/8"	22 1.5mm	07.2252.20					
07.0256.20	4.0mm	20.0cm/8"	22 2.0mm	07.2256.20					
07.0261.25	2.5mm	25cm/10"	26 1.2mm						
07.0264.30	3.0mm	30cm/12"	30 1.5mm	07.2264.30					
07.0268.30	4.0mm	30cm/12"	30 2.0mm	07.2268.30					

L = Working Length

07.2290.30

07-005

www.hilbro.net

07

EN Infiltration & Lipo-Suction Cannulas
DE Infiltration & Lipo-Suction Cannulas
ES & de la infiltración; Cánulas de Lipo-Suction

LIPO-SUCTION CANNULAS



22 holes,
Ø 1.5mm each

For Infiltration,
10 holes,
Ø 1.0mm each

Handle conn.	Ø	(L)	Luer-Lock	Handle conn.	Ø	(L)
07.0300.20	3.0mm	20cm/8"	07.2300.20	07.0310.20	2.0mm	25cm/10"
				07.0314.30	3.0mm	25cm/10"

07.0314.25

07-006

www.hilbro.net

L = Working Length

Handles & Adaptor for Cannulas **EN**
 Behandelt & Adapter für Cannulas **DE**
 Maneja el & Adaptador para las cánulas **ES**

07

LIPO-SUCTION CANNULAS



Ø L

07.0324.13 13.0mm 13cm/5¼"



With Hole for interrupting the suction process

Ø L

07.0334.13 13.0mm 13cm/5¼"

Adaptor Luer-Lock



L

07.2344.05 5.0cm/2"

L = Working Length

07.2344.05

07-007

www.hilbro.net

07

- EN Lipo-Suction Cannulas with Thumb Plate
- DE Fettabsaugung Cannulas mit Daumen-Platte
- ES Cánulas del Lipo-Suction con la placa del pulgar



LIPO-SUCTION CANNULAS



One central hole



Two central hole



Ø 13mm **CASPAROTTI**



Ø 13mm **CASPAROTTI**

	Ø	L		Ø	L
07.0354.15	3.0mm	15.0cm/6"		07.0380.15	4.0mm 15.0cm/6"
07.0358.15	6.0mm	15.0cm/6"		07.0380.25	4.0mm 25cm/10"
07.0362.25	4.0mm	25cm/10"		07.0384.30	6.0mm 30cm/12"
07.0366.25	6.0mm	25cm/10"			
07.0370.30	9.0mm	30cm/12"			


07.0384.30

07-008

www.hilbro.net

L = Working Length

LIPO-SUCTION CANNULAS



Ø 13mm **CASPAROTTI**

Ø 13mm **CASPAROTTI**

	Ø	L
07.0394.15	3.0mm	15cm/6"
07.0404.15	6.0mm	15.0cm/6"
07.0404.25	6.0mm	25cm/10"

L = Working Length

07.0404.25

07-009

www.hilbro.net